


Constant monitoring of the condition of patients in intensive care units






1: 1 h


Pleth 12.0 Pleth 100 98 PR 60 PR 60



Continuous recording of all sources 24/7

Archived video records are stored for 30 days for each connected video source







Patient history

Quick access to any event



Intuitive player

Synchronous playback of sources in any configuration









Archived post: 21 Nov 2020 11:23




Observation post


Doctors on duty have the opportunity to monitor changes and provide the necessary medical assistance in a timely manner













Teleconsultations and communication with relatives

Emergency button and all the benefits of telemedicine in one platform



Integration with information systems

Data exchange with the intensive care and anesthesiological

Mm Nº1 3 4. 5.

№ №2 2 2.2. №3 (Свободна) №4 (Свободна)

⊭ №1 Иванов Д. К

№2 (Свободна)

🛤 №1 Иванова 3. Я

I= №3 Тестин И. Д.

та 7 Иванов Д. К.


🚐 9 Романов П. А.


№ №1 Чапай В. И.

ОБЩЕЕ

information system







Video analytics

Development of a system for recognizing events and actions of personnel and patients

Solution for ICU

